



**ITABLOC®**



# **Boccaglio per gastroenterologia**

BITEBLOCK FOR GASTROENTEROLOGY

**Meditalia®**



**Meditalia**®

**BOCCAGLI PER GASTROENTEROLOGIA**  
BITEBLOCKS FOR GASTROENTEROLOGY

**Cintino premontato ad elasticità differenziata**  
*Pre-assembled headstrap with differentiated elasticity*

**Struttura rigida**  
*Sturdy construction*



**Abbassalingua munito di grip**  
*Tongue depressor with grip*

**Ampi accessi laterali**  
*Large side access*

**Raccordo in soft laprene atraumatico**  
*Nasal connector in soft and atraumatic laprene*

**Stampo cromato garantisce angolature smusse**  
*Chromed mold guarantees smooth angles*

**Products info**

References	Oral diameter	Material	Colour	Headstrap
MED/102	21x25 mm	polypropylene	yellow	•
MED/102/B	21x25 mm	polypropylene	yellow	•
MED/102/A	21x25 mm	polypropylene 	transparent	•
MED/102/A/B	21x25 mm	polypropylene 	transparent	•
MED/106	18x20 mm	polypropylene	purple	•
MED/106/B	18x20 mm	polypropylene	purple	•
MED/106/P	16x20 mm	polypropylene	orange	•
MED/106/P/A	16x20 mm	polypropylene	transparent	•